




# Endoscopes for Veterinary Laparoscopy

## Laparoscopes



- ★ 5mm diameter - 330mm length
- ★ 10mm diameter - 330mm length



### Trocars and cannulas

- ★ for 5mm endoscopes
- ★ for 10mm endoscopes



### Veress needles

- ★ 2mm Veress needles 10mm length with luer lock and stylet for insufflation

For **small animals** use the best scope is 5mm - 0 degree view. It is suitable for most of laparoscopy procedures.

For **large animals** use the best scope is 10mm - 0 degree or 30 degree view. To use it with you will need trocar and cannula. Veress needles is optional.

For best image results we recommend to use XENON light sources.

## Instruments and accessories

### Scissors



### Forceps

#### Biopsy forceps with spike



#### Punch biopsy forceps




#### Babcock grasping forceps



Other types of forceps available

### CO2 INSUFFLATORS



CO2 insufflator is essential part of laparoscopic equipment set. With this device you will be able to expand cavity and visualize it.

Please contact us for more information:

**MEDIT ENDOSCOPES, TEL (403) 242-0363  
FAX (613) 822-5255, www.endoscopy4you.com**