


# Universal Multi - Purpose Endoscope




2.7mm diameter - 188mm length, 30 degree view

## Why multi purpose?

The Multi-Purpose Rigid Endoscope is well suited to a broad range of procedures including: Rhinoscopy, Cystoscopy in female dogs, Urethroscopy, Otoscopy, Arthroscopy in dogs, Vaginoscopy, Laparoscopy, Thoracoscopy. The 2.7mm Multi-Purpose Rigid Endoscope is also ideal for Avian endoscopy.


## Operational sheath with 5 fr working channel

Operational sheath with working channel and two stopcocks 11FR and 14.5FR



Main use: Cystoscopy in female, Ectopic ureters diagnosis, Sinuscopy, Rhinoscopy with instruments use.

Optional: blunt obturator



## Flexible instruments





- ★ Foreign body forceps
- ★ Biopsy forceps



## Arthroscopy sheath (multi - purpose sheath)

This sheath is used for arthroscopy, sinuscopy, rhinoscopy, vaginoscopy, cystoscopy, thoracoscopy, otoscopy, etc.



Multi-purpose endoscope set is a good start for any small animal or exotic practitioner. Your diagnostic abilities will greatly increase with this affordable investment.