

Visual Systems for Veterinary Endoscopy

VIDEO CAMERAS

STAND ALONE VIDEO CAMERAS



MK 1 - low cost video camera system with high quality videoimage. This camera offer over 450 lines of resolution and have a most common capabilities. MK series cameras has superior image then other 1 chip camera brands on the market for medical or veterinary endoscopy.



MK 3 - high end camera system that give you a great opportunity to use endoscopy to the maximum extent! This camera give you the same good image then any other MK series cameras BUT it also give you CAPTURE and PROCESSING images capabilities. You will be able to save captured images on MEDIA card and use it it marketing your services to the clients and for hystory or later analysis.

COMBY VIDEO AND LIGHT SOURCE SYSTEM



OPTIC MASTER - comby low cost camera and l light source system. Include video camera with 410 lines of resolution and halogen light source. This is good economy system is great for low budget practitioners . This camera will work great for : Otoscopy, Sinuscopy, Urethrascopy Oral cavity scopy, vaginal scopy , rectoscopy. It can be used for arthroscopy and laparoscopy but we recomend MK systems for this applications.

Options:

- ★ Light source 250 watts, camera and freezing board (freezing images in memory)
 - ★ Ligth source 150 watts, camera without freezing board (most affordable option)
-
- ★ Hangun type camera head
 - ★ Straight type camera head

Visual Systems for Veterinary Endoscopy

MONITORS



MEGA MONITORS

Our monitors give 480 lines of resolution in a compact disain.

- ★ 14 inch
- ★ 15 inch

Monitor is essential part of you visual system. There is no excuse to use TV instead of a professional monitor. Dont save on monitors - save on something else! (the only exeption is to use TV with a optic master video camera)

LIGHT SOURCES



XENON LIGHT SOURCE - XE 300

- ★ Bright, white, intense light of 300 watts for ALL endoscopy applications.
- ★ Xenon light highly recomended for laparoscopy and large animal arthroscopy.



HALOGEN LIGHT SOURCE - CLH

- ★ Halogen, yellow, bright light
- ★ Good for otoscopy, rhinoscopy, cystoscopy, arthroscopy, etc. For laparoscopy and arthroscopy 250 watts light recomended.

MODELS:

- ★ CLH150 - 150 watts light
- ★ CLH250 - 250 watts light

ACCESSORIES

CAMERA COUPLER



LIGHT CABLE

Please contact us for more information:

MEDIT ENDOSCOPES, TEL (403) 242-0363
FAX (613) 822-5255, www.endoscopy4you.com